

Diffused Photon NIR Tomography for Traumatic Brain Injury

Portable FD-fNIRs system concept for early TBI screening

Thomas Molnar | Ryan Vo | Nathan Huynh

Group 8 | ENGR113 | Prof. Daryoush

Introduction

2

3,800,000


TBIs each year in the US

10%

Due to sports related injuries


492,000+

U.S. service members
diagnosed
with TBI from 2000-2023






- Magnetic Resonance Imaging (MRI)
- Computed Tomography (CAT Scan)
- X-Ray



Significance of Hemoglobin



- Protein that delivers oxygen
- Absorption spectra
- In healthy brains
 - About even oxygenated to deoxygenated
- Impact damages small blood vessels
 - Oxygenated hemoglobin leaks out
 - Abnormal ratio


Objectives

5

- Portable
- Fast
- Noninvasive
- Commercially available
- Implement
 - Transmitter
 - Receiver
- Calculations
- Imaging




Why Frequency-Domain Functional Near-Infrared Spectroscopy⁶



- Continuous Wave
 - Constant intensity
 - Measures change in intensity
- Time Domain
 - Short light pulses
 - Measures arrival time distributions
- Frequency Domain
 - Modulates emitted light intensity
 - Measures both intensity and phase

Why Three Wavelengths?


7



- Absorption of oxygenated and deoxygenated hemoglobin
 - Deoxygenated
 - 600-800nm
 - Isosbestic point
 - Approximately 800nm
 - Absorbed identically
 - Oxygenated
 - 850 nm
- Modified Beer-Lambert Law

Approach for the Project

8



FD-fNIR is a screening tool

Provides early information about possible brain injury and oxygenation changes; not intended to replace CT or MRI.

Accuracy depends on signal quality


Sensor placement, hair, motion, and poor skin contact can affect the received signal.

Source-detector spacing is a tradeoff

Greater spacing improves penetration depth, but reduces received signal strength.




General System Block Diagram




Transmitter Subsystem

10



Receiver Subsystem

11



Photodetector: Photodiodes + TIA + RF Switch + Splitter

Phase/Amplitude Detector: RF Amp, detectors, REF frequency

Data Acquisition: ADC + ESP32

Transmitter Subsystem Formula Map

RF Modulation

$$\Delta P_{\text{opt},\lambda} = \eta_{L,\lambda} I_{\text{RF},\lambda}$$

$$\omega = 2\pi f$$

$$I_{\text{laser},\lambda}(t) = I_{\text{bias},\lambda} + I_{\text{RF},\lambda} \cos(\omega t)$$

Sets the modulation frequency and RF current swing.

DC Bias / Linearity

$$\eta = \frac{P_{\text{optical}}}{P_{\text{elec}}} \times 100\%$$

$$I_{\text{bias},\lambda} - I_{\text{RF},\lambda} > I_{\text{th},\lambda}$$

$$P_{\text{elec},\lambda} = V_{F,\lambda} I_{F,\lambda}$$

Checks that the laser stays above threshold and in the linear modulation region.

Laser Output Stage

$$P_{\text{dBm}} = 10 \log_{10}(P_{\text{mW}})$$

$$\Delta P_{\text{opt},\lambda} = m_{\lambda} P_{\text{DC},\lambda}$$

$$P_{\text{tx},\lambda} = P_{\text{DC},\lambda} + \Delta P_{\text{opt},\lambda} \cos(\omega t)$$

Relates modulation depth to the optical power swing at 670, 800, and 850 nm.

Wavelength	Modulation Depth	Current Swing	Limit Check	Avg. dBm	Ptx (t) at 500 MHz
670nm	20%	4.50–6.00 mA	3.5<I<7 mA	0 dBm	1.5+0.30cos(ω t) mW
800nm	20%	4.50–5.50 mA	4<I<8 mA	–3.01 dBm	0.5+0.10cos(ω t) mW
850nm	50%	4.00–9.00 mA	1.5<I<12 mA	3.01 dBm	2.0+1.00cos(ω t) mW

$\omega=3.142\times10^9$ rad/s - 500MHz

Wavelength	VF	IF	Pelec =VF IF	Optical efficiency
670nm	2.1 V	4 mA	8.4mW	11.9%
800nm	1.9 V	5.0 mA	9.5 mW	5.3%
850nm	1.8 V	6.5 mA	11.7 mW	17.1%

Transmitter Calculations Optical Launch to Tissue

Optical Source Intensity

$$I_{0,\lambda} = \frac{P_{\text{tx},\lambda}}{A_{\text{beam},\lambda}}$$

Baseline optical intensity launched into tissue.

Effective Attenuation

$$\mu_{\text{eff},\lambda} = \sqrt{3\mu_{a,\lambda}(\mu_{a,\lambda} + \mu'_{s,\lambda})}$$

Combines absorption and scattering into intensity loss.

Tissue Propagation

$$I_{\lambda} = I_{0,\lambda} e^{-\mu_{\text{eff},\lambda} d}$$

Intensity decreases through forehead and brain tissue.

Receiver Input

$$P_{\text{det},\lambda} \propto I_{\lambda}$$

The detector input follows the tissue-reduced intensity.

Receiver Subsystem Signal Chain & Formula Map

Updated transmitter input

$$P_{tx,\lambda}(t) = P_{DC,\lambda} + \Delta P_{opt,\lambda} \cos(\omega t)$$

Detected power split

$$P_{det,DC,\lambda} = P_{DC,\lambda} e^{-\mu_{eff,\lambda} d}$$

$$\Delta P_{det,\lambda} = \Delta P_{opt,\lambda} e^{-\mu_{eff,\lambda} d}$$

DC path

$$P_{DC,\lambda} \rightarrow P_{det,DC,\lambda} \rightarrow I_{avg,\lambda} \rightarrow P_{noise,\lambda}$$

$$I_{avg,\lambda} = R_{D,\lambda} P_{det,DC,\lambda}$$

Used for average photodiode current and receiver noise.

AC modulation path

$$\Delta P_{opt,\lambda} \rightarrow \Delta P_{det,\lambda} \rightarrow V_{AC,\lambda} \rightarrow P_{sig,\lambda} \rightarrow SNR_\lambda$$

$$I_{AC,\lambda} = R_{D,\lambda} \Delta P_{det,\lambda}$$

Used for recovered modulation signal, IQ amplitude, signal power, and SNR.

Receiver Calculations Noise Power & SNR

Noise Current Terms

$$i_{\text{shot},\lambda}^2 = 2qI_{\text{avg},\lambda}B$$

$$i_{\text{RIN},\lambda}^2 = \text{RIN} \cdot I_{\text{avg},\lambda}^2 B$$

$$P_{\text{thermal}} = k_B T B$$

Noise grows with current, bandwidth, and electronics.

Total Receiver Noise

$$i_{\text{noise},\lambda} = \sqrt{i_{\text{shot}}^2 + i_{\text{RIN}}^2 + i_{\text{thermal}}^2}$$

$$P_{\text{noise},\lambda} = \frac{(i_{\text{noise},\lambda} Z_T)^2}{2R_L}$$

Converts total current noise into receiver noise power.

Signal Reliability

$$\text{SNR}_\lambda = \frac{P_{\text{sig},\lambda}}{P_{\text{noise},\lambda}}$$

$$P_{\text{sig},\lambda} > P_{\text{noise},\lambda}$$

Main metric: signal reliability and clear recovery.

IQ Demodulation & Signal Recovery

Received RF Signal

$$s_{\lambda}(t) = A_{\lambda} \cos(\omega t + \phi_{\lambda})$$

$$\omega = 2\pi f$$

Detected optical signal keeps amplitude and phase after tissue propagation.

I/Q Recovery

$$I_{\lambda} = \text{LPF}\{s_{\lambda}(t)\cos(\omega t)\}$$

$$Q_{\lambda} = \text{LPF}\{s_{\lambda}(t)\sin(\omega t)\}$$

$$A_{\lambda} = \sqrt{I_{\lambda}^2 + Q_{\lambda}^2}$$

$$\phi_{\lambda} = \tan^{-1}\left(\frac{Q_{\lambda}}{I_{\lambda}}\right)$$

Demodulation extracts baseband components and reconstructs amplitude and phase.

How to Use Them

$$I_{\text{rec}, \lambda} \propto A_{\lambda}$$

Amplitude is used for intensity and MBLL. Phase is handled separately for delay and path information.

Calculations Intensity-Based MBLL

Optical Density from Intensity

$$\Delta OD_{\lambda} = \log_{10} \left(\frac{I_{0,\lambda}}{I_{\lambda}} \right)$$

I_0 is baseline intensity; I_{λ} is measured received intensity.

Three-Wavelength Measurement

$$\lambda = 670 \text{ nm}, 800 \text{ nm}, 850 \text{ nm}$$

Each wavelength gives a separate optical-density equation.

Modified Beer–Lambert Law

$$\Delta OD_{\lambda} = \epsilon_{\text{HbO}_2, \lambda} \Delta[\text{HbO}_2] L_{\text{eff}, \lambda} + \epsilon_{\text{Hb}, \lambda} \Delta[\text{Hb}] L_{\text{eff}, \lambda}$$

MBLL is based on intensity change, not phase directly.

Matrix / Output

$$\Delta OD = E \Delta C L_{\text{eff}}$$

$$\Delta[\text{HbO}_2], \Delta[\text{Hb}]$$

Outputs: oxygenated and deoxygenated hemoglobin changes.

Calculations & Diagram

Wavelength Set

$$\lambda = 670, 800, 850 \text{ nm}$$

Selected wavelengths sample different hemoglobin responses.

Intensity $\rightarrow$ Optical Density


$$\Delta OD_{\lambda} = \log_{10} \left(\frac{I_{0,\lambda}}{I_{\lambda}} \right)$$

Each wavelength produces an optical-density change.

MBLL Estimation

$$\Delta OD_{\lambda} \rightarrow \Delta[\text{HbO}_2], \Delta[\text{Hb}]$$

Relative intensity changes across wavelengths support screening interpretation.



Hb and HbO₂ response curves with 670, 800, and 850 nm marked.

Core equations

$$\mu_{\text{eff},\lambda} = \sqrt{3\mu_{a,\lambda}(\mu_{a,\lambda} + \mu'_{s,\lambda})}$$

$$L_{\text{tissue,dB},\lambda} = 4.343 \mu_{\text{eff},\lambda} d$$

$$\Delta OD_{\lambda,d} = \log_{10} \left(\frac{I_{0,\lambda,d}}{I_{\lambda,d}} \right)$$

$$\uparrow \mu_a \Rightarrow \uparrow \mu_{\text{eff}} \Rightarrow \uparrow L_{\text{tissue}} \Rightarrow \downarrow I_{\lambda}$$

d1 / d2 table for μ_a and $\mu's$

Path	d	μ_a estimate	$\mu's$ estimate	Use
d1	2 cm	from $\Delta OD_{\lambda,d1}$	article/model fit	shallow clue
d2	4 cm	from $\Delta OD_{\lambda,d2}$	article/model fit	deeper clue

Loss for two depths

λ	μ_{eff}	L @ d1=2 cm	L @ d2=4 cm	Gain
670	2.31	20.1 dB	40.1 dB	+20.1 dB
800	2.04	17.7 dB	35.4 dB	+17.7 dB
850	2.17	18.8 dB	37.7 dB	+18.8 dB

Gain = less loss using the shallow 2 cm path instead of the deeper 4 cm path.

Leakage location clue

- Blood leakage can increase local absorption μ_a .
- Higher μ_a increases μ_{eff} and lowers detected intensity.
- Large ΔOD at d1 suggests a shallower absorption change.
- Large ΔOD at d2 suggests a deeper path is affected.
- Use 670, 800, 850 nm together to check hemoglobin response.

Formula used

$$P_{det,DC,\lambda} = P_{DC,\lambda} e^{-\mu_{eff,\lambda} d}$$

$$\Delta P_{det,\lambda} = \Delta P_{opt,\lambda} e^{-\mu_{eff,\lambda} d}$$

The DC part sets average receiver level; the AC part sets recovered modulation signal.

d = 2 cm

λ	Pdet,DC	ΔP_{det}
670	14.779 μW	2.956 μW
800	8.454 μW	1.691 μW
850	26.073 μW	13.037 μW

Meaning

- d = 2 cm gives stronger receiver input.
- d = 4 cm probes deeper but has much more attenuation.
- ΔP_{det} is the modulated signal that supports IQ and SNR.

d = 4 cm

λ	Pdet,DC	ΔP_{det}
670	0.146 μW	0.029 μW
800	0.143 μW	0.029 μW
850	0.340 μW	0.170 μW

Receiver equations

$I_{avg,\lambda} = R_D P_{det,DC,\lambda}$

$I_{AC,\lambda} = R_D \Delta P_{det,\lambda}$

$V_{out,DC} = I_{avg} Z_T, \quad V_{AC} = I_{AC} Z_T$

$P_{sig,\lambda} = \frac{V_{AC,\lambda}^2}{2R_L}$

Constants: RD = 0.65 A/W, ZT = 794 Ω, RL = 50 Ω

d = 2 cm

λ	I_{avg}	$V_{out,DC}$	I_{AC}	V_{AC}	P_{sig}
670	9.606 μA	7.628 mV	1.921 μA	1.526 mV	-46.33 dBm
800	5.495 μA	4.363 mV	1.099 μA	0.873 mV	-51.18 dBm
850	16.947 μA	13.456 mV	8.474 μA	6.728 mV	-33.44 dBm

d = 4 cm

λ	I_{avg}	$V_{out,DC}$	I_{AC}	V_{AC}	P_{sig}
670	0.095 μA	75.15 μV	0.019 μA	15.03 μV	-86.46 dBm
800	0.093 μA	73.77 μV	0.019 μA	14.75 μV	-86.62 dBm
850	0.221 μA	175.42 μV	0.110 μA	87.71 μV	-71.14 dBm

Optical density from intensity

$$\Delta OD_{\lambda} = \log_{10} \left(\frac{I_{0,\lambda}}{I_{\lambda}} \right)$$

Modified Beer-Lambert Law

$$\Delta OD_{\lambda} = \epsilon_{HbO_2,\lambda} \Delta[HbO_2] L_{eff,\lambda} + \epsilon_{Hb,\lambda} \Delta[Hb] L_{eff,\lambda}$$

Matrix form

$$\Delta OD = E \cdot \Delta C \cdot L_{eff}$$

Extinction coefficient table

λ	ϵ_{HbO_2}	ϵ_{Hb}
670	0.000294	0.002795
800	0.000816	0.000762
850	0.001058	0.000691

Expected output

- $\Delta[HbO_2]$
- $\Delta[Hb]$
- Oxygenation trend from intensity changes.

SNR equation

$$SNR_{dB} = P_{sig, dBm} - P_{noise, dBm}$$


The signal power uses the AC recovered modulation signal from ΔP_{det} .

d = 2 cm				
λ	Pnoise	Psig	SNR	
670	-135.95 dBm	-46.33 dBm	89.62 dB	
800	-138.00 dBm	-51.18 dBm	86.82 dB	
850	-133.56 dBm	-33.44 dBm	100.12 dB	

- Meaning
- Higher SNR supports clearer intensity recovery.
 - 2 cm has higher SNR because less tissue loss occurs.
 - 4 cm has lower SNR but gives deeper-path information.

d = 4 cm				
λ	Pnoise	Psig	SNR	
670	-143.63 dBm	-86.46 dBm	57.17 dB	
800	-143.64 dBm	-86.62 dBm	57.02 dB	
850	-143.38 dBm	-71.14 dBm	72.24 dB	

- Portable FD-fNIR system designed for early TBI screening.
- Intended as a real-time screening tool, not a replacement for CT or MRI.



Thank You

Questions?

Diffused Photon NIR Tomography
for Traumatic Brain Injury

Thomas Molnar | Ryan Vo | Nathan Huynh

ENGR113 | Prof. Daryoush

Detailed Calculations - Transmitter Frequency & Optical Power

Modulation frequency

$$\omega = 2\pi f = 2\pi(500 \times 10^6) = 3.142 \times 10^9 \text{ rad/s}$$

Calculated transmitter optical power

λ	Popt	Popt dBm
670 nm	5.0 mW	6.99 dBm
800 nm	1.65 mW	2.17 dBm
850 nm	1.85 mW	2.67 dBm

Optical output power conversion

$$P_{\text{dBm}} = 10\log_{10}(P_{\text{mW}})$$

$$670 \text{ nm} : 10\log_{10}(5.0) = 6.99 \text{ dBm}$$

$$800 \text{ nm} : 10\log_{10}(1.65) = 2.17 \text{ dBm}$$

$$850 \text{ nm} : 10\log_{10}(1.85) = 2.67 \text{ dBm}$$

Meaning

- dBm conversion helps compare optical power with dB losses.
- This supports the transmitter-to-receiver link budget.

Detailed Calculations - RF Current & Laser Linearity

RF modulation current

$$I_{\text{RF}} = \sqrt{\frac{2P_{\text{RF}}}{R_F}}, \quad P_{\text{RF}} = 1.0 \text{ mW}$$

$$R_F = 5\Omega: I_{\text{RF}} = \sqrt{\frac{2(1.0 \times 10^{-3})}{5}} = 20.0 \text{ mA}$$

$$R_F = 10\Omega: I_{\text{RF}} = \sqrt{\frac{2(1.0 \times 10^{-3})}{10}} = 14.1 \text{ mA}$$

RF current range: 14.1–20.0 mA

Laser current and linearity

$$I_{\text{laser}, \lambda} = I_{\text{bias}, \lambda} + I_{\text{RF}, \lambda} \cos(\omega t)$$

$$I_{\text{bias}, \lambda} - I_{\text{RF}, \lambda} > I_{\text{th}, \lambda}$$

Numbers to check

- $P_{\text{RF}} = 1.0 \text{ mW}$
- $R_F = 5\text{--}10 \Omega$
- $I_{\text{RF}} = 14.1\text{--}20.0 \text{ mA}$
- I_{bias} and I_{th} come from the final laser datasheet

Meaning

- RF current controls the strength of intensity modulation.
- Bias condition prevents clipping and threshold drop-out.

Detailed Calculations - Transmitter Electrical Power & Efficiency

Electrical power

$$P_{elec} = V_F I_F$$

- Use final datasheet values for V_F and I_F .
- Expected electrical power stays in the mW range.

Laser efficiency

$$\eta = \frac{P_{optical}}{P_{elec}} \times 100\%$$

$$670 \text{ nm} : \eta = \frac{5.0}{25.0} \times 100\% = 20.0\%$$

$$800 \text{ nm} : \eta = \frac{1.65}{20.0} \times 100\% = 8.25\%$$

$$800 \text{ nm} : \eta = \frac{1.65}{13.6} \times 100\% = 12.1\%$$

$$850 \text{ nm} : \eta = \frac{1.85}{17.6} \times 100\% = 10.5\%$$

$$850 \text{ nm} : \eta = \frac{1.85}{13.6} \times 100\% = 13.6\%$$

Expected values

λ	P_{opt}	P_{elec}	η
670	5.0 mW	25.0 mW	20.0%
800	1.65 mW	13.6-20.0 mW	8.25-12.1%
850	1.85 mW	13.6-17.6 mW	10.5-13.6%

Meaning

- Electrical power estimates laser drive requirements.
- Efficiency shows how much electrical input becomes optical output.

Detailed Calculations - Optical Modulation & Transmitted Power

Optical modulation depth

$$\Delta P_{\text{opt},\lambda} = \eta_{L,\lambda} I_{\text{RF},\lambda}$$

- Plug in $I_{\text{RF}} = 14.1\text{--}20.0\text{ mA}$.
- Plug in η_L from each selected laser diode datasheet.

Transmitted optical waveform

$$P_{\text{tx},\lambda} = P_{\text{DC},\lambda} + \Delta P_{\text{opt},\lambda} \cos(\omega t)$$

- $P_{\text{DC},670} = 5.0\text{ mW}$
- $P_{\text{DC},800} = 1.65\text{ mW}$
- $P_{\text{DC},850} = 1.85\text{ mW}$

Expected output

- A DC optical output plus RF optical modulation.
- This modulated light enters the forehead/brain tissue.
- Exact ΔP_{opt} needs final slope efficiency.

Meaning

- This connects transmitter hardware to the optical signal.
- The receiver later detects the attenuated version of this signal.